

<b>Prepared by:</b> Steve Brown	<b>SOP – TISSUE CULTURE WORKUP</b>
<b>Revised by:</b>	<b>April 25, 2006</b>

**PURPOSE:**

To establish a standard operating procedure for tissue culture infection and measurement of tau.

**PROCEDURE:**

**VIRUS PREPARATION**

1. Split 293T cells to 10cm dishes so that they are semi-confluent the day of transfection.
  - 1.1. Usually a confluent plate split 1:6 the day before transfection is about correct.
  - 1.2. Plates should contain 8ml DMEM/10%FCS/1%pennicillin-streptomycin-glutamine
2. Make DNA mix – plasmids, 100ul CaCl<sub>2</sub> SOLUTION, H<sub>2</sub>O to 500ul.
  - 2.1. CaCl<sub>2</sub> SOLUTION – 2M CaCl<sub>2</sub>, tissue culture grade, filter-sterilized.
  - 2.2. Plasmids – 14ug coding plasmid, 10ug p8.91 packaging plasmid, 6ug pMD2G coat plasmid.
    - 2.2.1. Plasmids should be of high purity: Quiagen maxi-prep, JetStar maxi-prep, or CsCl<sub>2</sub>-purified, dissolved in H<sub>2</sub>O.
  - 2.3. Water used for transfection should be buffered with 2mM HEPES pH 7.05.
3. In 4-10ml polypropylene tube, place 0.5ml 2xHeBS.
  - 3.1. 2xHeBS – 6mM dextrose, 137mM NaCl, 5mM KCl, 0.75mM Na<sub>2</sub>HPO<sub>4</sub>\*7H<sub>2</sub>O, 25mM HEPES (free acid). Then adjust to pH 7.05 with NaOH
  - 3.2. pH of HeBS is absolutely critical. 7.05 is ideal, less than 7.0 or more than 7.1 is useless.
4. Add DNA to HeBS solution in fashion specified.
  - 4.1. All solutions should be at room temperature.
  - 4.2. Bubble air into HeBS with pipet. Add DNA dropwise with pipet tip, letting it slide down side of bubbling pipet.
  - 4.3. Vortex 10 seconds.
  - 4.4. Incubate exactly 30 minutes at room temperature.
5. A fine white translucent precipitate should form. Add dropwise to medium.
  - 5.1. Separate mix per plate; do not do in batch.
  - 5.2. Mix precipitate and medium by rocking plate back and forth, as for spreading cells.
6. After 24 hours of incubation in standard incubator, change medium.
  - 6.1. Rinse cells 1x with PBS.
  - 6.2. Replace medium with 10ml fresh medium.
7. After 24 hours, harvest medium containing virus, and replace with 10ml fresh medium.
  - 7.1. Store viral supernatant on ice.
8. After another 24 hours, harvest supernatant again and pool with first medium.
  - 8.1. Plates can now be discarded.
9. Spin supernatant 2k 5 minutes to pellet debris.
10. Filter through 0.2-0.4micron filters. Store in aliquots at -80.
11. Supernatant can be concentrated by ultracentrifugation, 30,000 rpm 90 minutes in mini-ultracentrifuge or comparable.
12. Resuspend overnight in desired volume of medium by medium-speed vortexing in automatic vortexer.
  - 12.1. 10x concentration is ideal. More than this is toxic.

**INFECTION AND MEASUREMENT**

1. Thaw primary fibroblasts.
  - 1.1. Unfreeze quickly in 37 degree waterbath, and dilute into 10ml PBS.
  - 1.2. Spin 1000 rpm 5 minutes, and aspirate supernatant.
  - 1.3. Resuspend in 2ml growth medium, and plate to 4 cm. dish, if appropriate.
2. The day before infection, split confluent cells 1:4 to 2x4cm dishes.
3. The day of transfection, remove medium and add warmed viral supernatant.
  - 3.1. Don't forget to increase serum to 20% if appropriate.
  - 3.2. For immortalized cells, addition of polybrene to 1ug/ml final concentration is often helpful. To primary cells, rather add 8ug/ml protamine sulfate final concentration.
4. Change to fresh supernatant after 6-24 hours, and let cells recover overnight.
5. Repeat 1x-2x as desired. If cells become confluent during this time, split 1:2 again as necessary.
6. Rinse cells, and replace medium with fresh growth medium.
7. After 3 days, conduct measurement.
  - 7.1. Treat cells 15 minutes with INDUCTION MEDIUM.
    - 7.1.1. INDUCTION MEDIUM – growth medium + 100nM dexamethasone (from 1/10,000 stock in EtOH)
  - 7.2. Rinse cells 2x with 1ml PBS.
  - 7.3. Replace medium with 3ml COUNTING MEDIUM.
    - 7.3.1. COUNTING MEDIUM – DMEM high-glucose without phenol red, 20% FCS, 1% pen-strep-glu, 0.1mM luciferin (from 1000x stock in H<sub>2</sub>O).
  - 7.4. Place cells into machine.
8. Count in 5-minute bins for 5 days or until rhythms dampen to flatness.